


# Identification of Microprotein-Protein Interactions via APEX Tagging

*Qian Chu<sup>1</sup>, Annie Rathore<sup>1</sup>, Jolene K. Diedrich<sup>1,2</sup>, Cynthia J. Donaldson<sup>1</sup>, John R. Yates III<sup>2</sup>, Alan Saghatelian<sup>1</sup>*

<sup>1</sup>The Salk Institute for Biological studies, La Jolla, CA <sup>2</sup>The Scripps Research Institute, La Jolla, CA


## smORF/Microprotein



smORF	Length	Interacting protein	Biology
<b>Bacteria</b>			
SgrT	43 aa	Glucose transporter (PtsG)	Glucose metabolism
<b>Flies</b>			
Tal/Pri	11/32 aa	Ubr3	Development
<b>Mice</b>			
Mln	46 aa	Calcium transporter (SERCA1)	Muscle contraction and endurance
<b>Human</b>			
MRI/Cyren	69 aa	Ku70/Ku80	DNA repair
Nobody	68 aa	EDC4/Dcp1	mRNA decapping



- Microproteins are a new class of molecules that are encoded directly from small open reading frames (smORFs)
- Microproteins have fundamental roles in biology, including metabolism, apoptosis, and development
- Most microprotein functions rely on microprotein-protein interactions

## In situ APEX biotinylation




- An engineered ascorbate peroxidase (APEX) allows proximity biotin tagging in live cells
- APEX fusion to a microprotein of interest (MPOI) generates biotin phenoxyl radicals, and thus biotin labels surrounding proteins, which are likely microprotein interacting partners

## Case study I: MRI-APEX





- APEX fusion to MRI microprotein successfully biotinylates and enriches Ku70 and Ku80
- Compared to FLAG IP, most MS hits in MRI-APEX tagging are nuclear proteins, which is consistent with its nuclear localization
- In addition to the Ku70/Ku80 heterodimer, MRI-APEX tagging also reveals many other proteins in the NHEJ DNA repair pathway

## Case study II: C11orf98-APEX



- C11orf98-MP is an uncharacterized 123 aa microprotein encoded by C11orf98 smORF
- It has robust RNA-Seq and Ribo-Seq coverage, and is conserved between human and mouse, suggesting a potential function

Human	1	MGAPGGKINRPRT	ELKKLFKRRVLR	RRRLHR	VVGAVIDGGLI	TRHHLKRRSSAR	ANITLSGKKRR	70
Gorilla	1	MGAPGGKINRPRT	ELKKLFKRRVLR	RRRLHR	VVGAVIDGGLI	TRHHLKRRSSAR	ANITLSGKKRR	70
Chimpanzee	1	MGAPGGKINRPRT	ELKKLFKRRVLR	RRRLHR	VVGAVIDGGLI	TRHHLKRRSSAR	ANITLSGKKRR	70
Mouse	1	MAPPGGKINRPRT	ELKKLFKRRVLR	RRRLHR	VVGAVIDGGLI	TRHHLKRRSSAR	ANITLSGKKRR	70
Rat	1	MAPPGGKINRPRT	ELKKLFKRRVLR	RRRLHR	VVGAVIDGGLI	TRHHLKRRSSAR	ANITLSGKKRR	70



- APEX fusion to C11orf98 results in biotinylation and enrichment of many nuclear proteins
- Nucleolin (NCL) and nucleophosmin (NPM1) are the two most promising candidates as they have the highest fold enrichment and form a protein complex in the nucleus
- FLAG IP and immunofluorescence experiments validate the interaction between C11orf98-MP, NCL and NPM1

## Conclusion

- APEX tagging is able to biotinylate and enrich microprotein interacting proteins in live cells
- Compared to FLAG immunoprecipitation, APEX tagging demonstrates greater fold increase of bona fide interacting proteins, while decreasing background non-specific interactions
- The apparent improvement in the APEX data supports the application of APEX to the remaining uncharacterized microproteins

## Acknowledgements

Q.C. is a postdoctoral fellow funded by the George E. Hewitt Foundation for medical research. This study was supported by the NIH with an NCI Cancer Center Support Grant P30 (CA014195 MASS core, A.S.) and R01 (GM102491, A.S.), the Leona M. and Harry B. Helmsley Charitable Trust grant (#2012-PG-MED002, A.S.), and Dr. Frederick Paulsen Chair/Ferring Pharmaceuticals (A.S.).